

 labodam



info@labodam.com | www.labodam.com

204 L Flammable Storage Cabinet LD-LFSC-A12

204 L Flammable Storage Cabinet LD-LFSC-A12 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three-point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead-free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 204 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three-point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications

Model No.	LD-LFSC-A12
Capacity	204 L
Door Type	Dual Door Manual
Adjustable Shelf	1
Shelves Loading Capacity	100 kgs
Dimension (H x W x D)	1650 x 860 x 860 mm
Weight	202 kgs

labodam

Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA

United Kingdom

www.labodam.com // info@labodam.com