



*227 L Flammable Storage Cabinet*  
*LFSC-A13*

## Description :

227 L Flammable Storage Cabinet LFSC-A13 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

## Specifications :

Capacity	227 L
Door type	Dual door manual
Adjustable shelf	2
Shelves loading capacity	100 kgs
Dimension ( H x W x D )	1650 x 860 x 860 mm
Weight	173 kgs

## Features :

- Capacity - 227 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelf
- Adjustable spill retaining shelves

## Applications :

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)



**Labodam Equipment Ltd**

18a Melton Road Leicester LE4 5EA  
United Kingdom

[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)