



**83 L WEAK ACID AND ALKALI  
CABINET LD-LWAC-B12**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

## 83 L Weak Acid and Alkali Cabinet LD-LWAC-B12

83 L Weak Acid and Alkali Cabinet LD-LWAC-B12 is a single walled manual acid and alkali cabinet with 50 mm of insulating air space. Doors are reinforced with a three-point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers. Inner and outer epoxy resin coating provides durability and chemical resistance preventing corrosion and humidity.

### Features

- Capacity - 83 L
- Single walled manual
- 50 mm insulating air space
- Three-point latch and manual lock
- Three adjustable shelves

### Applications

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid, Carbonic acid, phosphoric acid, silicic acid, etc.

### Specifications

Model No.	LD-LWAC-B12
Capacity	83 L
Door Type	Single door manual
Adjustable Shelf	1
Dimension (H x W x D)	700 x 570 x 1815 mm
Weight	116 kgs



**Labodam Equipment Ltd**

18a Melton Road Leicester LE4 5EA  
United Kingdom

[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)