



*Variable Volume Multi Channel
Pipette VMP204L*

Description :

Specifications :

No of Channels	12
Volume Range	5-50 µL
Increment	0.5 µL
Test Volume	50 µL
	25 µL
	5 µL
Accuracy error	1.00 %
	1.50 %
	3.00 %
Precision error	0.50 %
	1.00 %
	2.00 %

Features :

Autoclavable at 121 °C

Material: Polycarbonate

Volume settings can be clearly read through digital display

Effortless pipetting accessibility with rotatable Dispensing head

Easy repair and maintenance with separate piston and tip cone assembly

Compatible with most universal tip

The logo for Labodam features the word "labodam" in a bold, black, sans-serif font. The letter "l" is significantly larger than the other letters. A green circular icon is positioned above the "a", consisting of a solid green circle with three small white dots above it, resembling a stylized atom or a drop.

Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA

United Kingdom

www.labodam.com // info@labodam.com