



*Variable Volume Multi Channel
Pipette VMP205L*

www.labodam.com // info@labodam.com

Description :

Specifications :

No of Channels	12
Volume Range	50-300 µL
Increment	5 µL
Test Volume	300 µL 150 µL 50 µL
Accuracy error	0.70 % 1.00 % 1.50 %
Precision error	0.25 % 0.50 % 0.80 %

Features :

Autoclavable at 121 °C

Material: Polycarbonate

Volume settings can be clearly read through digital display

Effortless pipetting accessibility with rotatable Dispensing head

Easy repair and maintenance with separate piston and tip cone assembly

Compatible with most universal tip



Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA

United Kingdom

www.labodam.com // info@labodam.com